South African Veterinary Laboratory Scientific Forum





Enquiry: Dr. J.G. Kangumba Ref.: 01/2011 Tel. No.: (018) 294 3232/3 Fax No.: (018) 293 2385 Email: <u>JKangumba@nwpg.gov.za</u>

OFFICE OF THE CHAIRPERSON

MEMORANDUM

- TO : SAVLSF DAFF-Approved Laboratories
- FROM : SAVLSF Chairperson
- DATE : 14 January 2011

SUBJECT : SADC HARMONIZED SOP FOR AI HA/HI (VERSION 2.0) AND THE USE OF 2 ANTIGENS

Some queries were received with regard to the use of 2 antigens.

This memo serves to clarify the correct interpretation of the relevant sections of the abovementioned SOP: Section 7 (7.1) bullet one and Section 11.

- 1. The use of 2 antigens shall mean:
 - a. $H5N_a$ and $H5N_b$ / $H7N_a$ and $H7N_b$ at the same time (i.e. a total of 4 antigens if testing sera for antibodies against both H5 and H7), or
 - b. 1 antigen (H5N_a / H7N_a) with subsequent testing (**compulsory**) of sera found positive on H5 / H7 with a <u>second</u> antigen H5N_b / H7N_b
- 2. The SOP has not been amended

Sincerely,

Dr. Joule Kangumba Chairperson: SAVLSF